

**Repeatability of granulation test**


Mesh	#1	#2	#3	#4
12 and up	5.6	0.8	1.2	1.3
12~16	0.5	0.9	1	1.2
16~24	0.6	0.8	1.2	1.4
24~32	0.7	0.8	0.9	1.1
32~42	1.6	1.7	1.9	1.8
42~60	5.9	4.3	4.8	3.5
60~80	9.6	8.5	8.5	6.6
80~115	13.2	15.6	13.4	12.8
115 and under	66.8	66.6	67	70.6
Average particle size*	135.6	135.7	138.3	136.9

\*Average particle diameter of the geometric mean

(Conditions)

Raw material	Sintered alumina (average particle size 40) 400g
Binder	5% PVA solution (#500) 25g
Inlet temperature	100°C
Binder liquid feed rate	12.4g/min
Binder spray times	6 times
Binder spray pressure	78kPa(0.8kg/cm <sup>2</sup> )
Nozzle height	25cm from microporous plate


The granulation process has many operation factors, the reproducibility depends on the skill level of the operation. The flow state of the granules has a significant impact on the test results. By adjusting the amount of hot air consistent flow conditions are achievable.

**Change of particle diameter**

The factors that influence the particle diameter are the binder liquid feed rate and the spray pressure, the former being the most influential. A higher binder amount will result in larger diameter particles.

(Conditions)

Raw material	Lactose(100 mesh under) 200g
Binder	70% Sorbitol solution
Inlet temperature	90°C
Binder liquid feed rate	12g/min
Binder spray times	6 times
Binder spray pressure	98kPa (1.0kg/cm <sup>2</sup> )
Nozzle height	25cm from microporous plate

**Repeatability of granulation test**

Particles generated by the pulvis mini bed are usually in the range of 0.1~1.5a. The particle size uniformity is lower than extrusion granulation and compression granulation methods.

The granularity consistency may be regulated by test conditions.

(Conditions)

Raw material	Lactose (100 mesh under) 200g
Binder	70% Sorbitol solution 7.3g
Inlet temperature	90°C
Binder liquid feed rate	12g/min
Binder spray times	7 times
Binder spray pressure	98kPa(1.0kg/cm <sup>2</sup> )
Nozzle height	22.5cm from microporous plate

**Example of implementation (Spray dryer ADL311SA)**

Sample name	Composition (%)	Inlet temp. (°C)	Outlet temp. (°C)	Dry air amount (m <sup>3</sup> /min)	Spray air pressure kPa(kg/cm <sup>2</sup> )	Sent amount of sample liquid (g/min)	Sample recovery rate (%)
Dextrin (solution)	10	150	80	0.4	98 (1.0)	6.1	66
Dextrin (emulsion)	40	150	80	0.4	98 (1.0)	5.1	63
Oxidized titanium (suspended liquid)	10	150	85	0.42	98 (1.0)	5.3	50
Soy sauce	50	130	75	0.36	98 (1.0)	5.1	60
Salt	10	145	85	0.38	98 (1.0)	5.3	52

**Repeatability of spray drying test (spray dryer ADL311SA)**

Test No.	Sample name	Sample density (%)	Drying conditions						Yield (g)	Recovery rate (%)	
			Inlet temp. (°C)	Outlet temp. (°C)	Dry air amount (m <sup>3</sup> /min)	Spray air pressure kPa(kg/cm <sup>2</sup> )	Test sample amount (g/min)	Sent amount of sample liquid (g/min)			
1	Coffee solution	5.00	150	75	0.45	147(1.5)	93.1	3.1	30	4.3	92.4
2	Coffee solution	5.00	150	75	0.45	147(1.5)	93	3.1	30	4	86
3	Coffee solution	5.00	150	75	0.45	147(1.5)	91.4	2.0	30	4	87.5
4	Coffee solution	5.00	150	75	0.45	147(1.5)	84.9	2.8	30	3.7	87.2
5	Coffee solution	5.00	150	75	0.45	147(1.5)	83.8	2.8	30	3.7	88.3

**Example of implementation (Pulvis mini spray GB-210A)**

Sample name	Sample density	Inlet temp.	Outlet temp.	Dry air amount	Spray air pressure	Sent amount of sample liquid	Recovery rate (%)
Dextrin (solution)	20% solution	140	85	0.48	147(1.5)	8.8	66
Drug suspension	10% suspension	145	80	0.42	196(2.0)	8.2	82
Black tea extract	20%solution	155	100	0.4	147(1.5)	7.8	72
Silica gel	20%solution	150	75	0.48	147(1.5)	12.6	70
Iron oxide	3%suspension	175	90	0.4	127(1.3)	9.5	75

**Example of implementation (Pulvis mini bed GB-210B)**

Sample	Binder			Test conditions							Results	
	Name	Weight (min)	Name	Concentration (%)	Spray amount (min)	Inlet temp. (°C)	Liquid sending rate (g/min)	Spray pressure kPa(kg/cm <sup>2</sup> )	Spray times (times)	Nozzle height (cm)	Average dia. (μm)	12~115 mesh recovery rate(%)
Silicon	200	PVA	5.0	77	125	15	59 (0.6)	4	27	339	58	
Oxidized iron	160	PVA	2.5	50	120	15	98 (1.0)	4	21	205	62	
Ceramics	200	PVA	3.0	106	120	15	78 (0.8)	3	22	404	82	
Alumina	160	PVA	3.0	60	110	15	59 (0.6)	4	22	311	88	
Silica	150	CMC	1.0	100	120	15	78 (0.8)	4	22	306	60	
Lactose	200	Sorbitol	70.0	10	100	14	98 (1.0)	4	25	390	80	
Black tea essence	250	Guar gum	0.5	24	85	6	59 (0.6)	10	28	333	77	
Grease containing powder	200	Glucose	30.0	11	85	4	59 (0.6)	7	22	236	82	

**Binder category and features**

Category	Features
Gelatin	Gelatin Low density and weak bonding strength. No need to heat.
Dextrin	While it has excellent disintegrating and formability, the binding strength is weak.
Potato starch	Good granulation properties and inexpensive. Used in the pharmaceutical and food sector.
Arsinic acid soda	Suitable as a binder for the high viscosity samples. Used primarily in the food sector.
Gum arabic	Warm and spray. Need large amount of binder.
CMC (Carboxymethyl cellulose)	High viscosity at low temperatures. High amount of powder remains.
HPC (hydroxypropyl cellulose)	Good cohesion and is suitable for hydrophilic material